



Specification Sheet

Revised 03/01/06

PRODUCT: New-Gar® deodorized 5000 Allicin powder garlic (the new & improved Nutra-Gar)		
PRODUCT CODE No.: NG5000P	Technical Administrator: Yang Xian	
PLANT SPECIES: <i>Allium Sativum</i>	<u>Yang Xian</u>	
Part used: Bulb	Complies with all current USP/NF	
SOURCE: China	Harvest Time : Sept/Oct	
Item	Specification	Method of Analysis

Appearance:

White to light tan powder

Visual

<u>Active Ingredients</u>	<u>Specification</u>	<u>Method of Analysis</u>
Allicin	Min. 5,000 ppm	C-18 HPLC
Alliin	Min 11,200 ppm	C-18 HPLC
Total Thiosulfinates	Min 5,000 ppm	C-18 HPLC
Gamma-Glutamylcysteines	Min 4,000 ppm	C-18 HPLC
Total Sulfur	VARIES	AOAC 990.28

Pesticide Evaluation

Organophosphorus Insecticides		GC-MS, USP 27, 22<561>, 2004
Diechlorvos	ND	
Fonofos	“	
Diazanon		
Parathion -methyl		
Chlorpyrifos-mythyl	“	
Pirimiphos-mythyl		
Malathion		
Parathion		
Chlorpyrifos	«	
Methidathion		
Ethion		
Carbophenothion		
Azinphos-methyl		
Phosalone		
Organochlorine & Pyrethroid Insecticides	ND	GC-MS, USP27, 22<561>, 2004
<i>alfar</i> -Hexachlorocyclohexane	“	
Hexachlorobenzene		
<i>bata</i> -Haxachlorocyclohexane		
Lindane		
<i>gamma</i> -Haxachlorocyclohexan		
Heptachlor		
Aldrin		
Heptachlor epoxide	ND	

Specification sheet

PRODUCT: Deodorized New-Gar ®5000 ppm allicin powder garlic (NG5000P)		
Item	Specification	Method of Analysis

<u>Pesticide Evaluation</u>	<u>Specification</u>	<u>Method of Analysis</u>
<i>o' p' -DDE</i>	ND	GC_USP27NF22/2004
alfar -Endosulfan		
Dieldrin		
<i>p' p' -DDE</i>		
Endrin		
Endosulfan		
<i>o' p' -DDT</i>		
Carbophenothion		
Organochlorine & Pyrethroid	“	
Insecticides		
<i>p' p' -DDT</i>		
Permethrin		
<i>Trans-Permethrin</i>		
Cypermethrin		
Fenvalerate		
Deltamathrin		
<u>Heavy metals</u>	<5mg/kg	USP method 231
Arsenic	<0.50mg/kg	
Cadmium	<0.50mg/kg	
Copper	<0.50mg/kg	
Lead	<0.50mg/kg	
Mercury	<0.50mg/kg	
Zinc	<0.50mg/kg	
<u>Bulk Density</u>	0.42-0.62 g/ml	
<u>Moisture</u>	7.0% max (although we aim for <6%)	A.D.O.G.A.
<u>IR Identification</u>	Spectrum matches reference IR	GIR001
<u>Granulation/Mesh</u>		
% on US #35	0% max	A.D.O.G.A.
% on US #45	1% max	A.D.O.G.A.
% on US # 80	15% max	A.D.O.G.A.
% thru US #100	50% max	A.D.O.G.A.
<u>Microbiological</u>		
Total plate Count	<10,000 cfu/g (Same as USP method)	AOAC 966.23
Yeast and Mold	<100 /g	FDA-BAM 7 th ed.
Coliforms	<20 /g	AOAC 966.23
E.Coli	negative	AOAC 966.23
Pseudomonas	negative	AOAC /BAM
Salmonella	negative	AOAC /BAM
Staph Aureus	negative	AOAC /BAM

Actives Testing : Each batch is tested for Allicin potential, *all other actives results are by input

Ingredients: 100% garlic, no chemicals or additives & NO ETO Treatment

pesticides & H/M are tested at US labs twice/yr : early & late crop

Storage Conditions & shelf life:

Store in a cool, dry place, do not freeze. Shelf life is 2-3 yrs depending on the value of allicin potential & based on our completed stability study.